



Bedside Cleaning Add Water Products



APPLI-KIT™ BEDSIDE PRECLEAN KIT

| Description | Code | |
|---|---|---|
| <p>For optimal rinsing and pre-cleaning of instruments at bedside</p> <p>Delivers Low Foaming D-Zyme™ enzymatic detergents containing 3 enzymes - Protease, Amylase and Lipase, PLUS Surfactant Detergents to break down proteinaceous bioburden during the preclean of surgical instruments</p> <ul style="list-style-type: none"> Green Friendly, 100% Medical Grade foam sponges, No Flaking or Scratching Biodegradable, non-toxic, odorless and non-irritating to the eyes, bacteriostatic <p>Stable wide base pouch for easy fluid transfer. Clear pouch bag for a direct view of the CONTENTS</p> <p>Add up to 500 ml of water to the POUCH then soak and squeeze the sponge 3 - 5 times to activate the D-Zyme™ enzymatic concentrate.</p> | |  <p>Appli-Kit™ Flat- Add Water to bag</p> |
| <p>Appli-Kit™ w/ D-Zyme™ Concentrate in stand-up pouch bag</p> <p>Options:</p> <ul style="list-style-type: none"> Appli-Kit™ Flat w/ D-Zyme™ Concentrate - Add Water to pouch (up to 500 ml) 100 ct / 1 case Appli-Kit™ Tubular w/ D-Zyme™ Concentrate - Add Water to pouch (up to 500 ml) 70 ct / 1 case | <p>MED-AFC500100</p> <p>MED-ATC500070</p> |  <p>Appli-Kit™ Tube - Add Water to bag</p> |



info@falconfireinc.com
www.falconfireinc.com

FALCONFIRE TECHNOLOGIES INC.



80 Bradford St., Suite WW23
Barrie, ON L4N 6S7

t: 866-635-2550



Bedside Cleaning Add Water Products



D-Zyme™ Concentrate Sponge

| Description | Code | |
|---|---|--|
| <p>For optimal rinsing and pre-cleaning of instruments at bedside</p> <p>Activate concentrate with up to 500 ml of water</p> <p>Delivers Low Foaming D-Zyme™ enzymatic detergents containing 3 enzymes - Protease, Amylase and Lipase, PLUS Surfactant Detergents to break down proteinaceous bioburden during the pre-cleaning of surgical instruments</p> <ul style="list-style-type: none"> Green Friendly, 100% Medical Grade foam, No Flaking or Scratching Biodegradable, non-toxic, odorless and non-irritating to the eyes, bacteriostatic <p>Options:</p> <ul style="list-style-type: none"> D-Zyme™ Concentrate Flat - Add to Water - 200 ct / 1 case D-Zyme™ Concentrate Tubular - Add to Water - 100 ct / 1 case | <p>MED-EWC805200</p> <p>MED-EWC804100</p> |  <p>D-Zyme™ Concentrate Flat</p>  <p>D-Zyme™ Concentrate Tubular</p> |

D-Zyme™ Concentrate with 500 ml Basin

| | | |
|--|---|--|
| <p>Add up to 500ml of water to the clear plastic basin, then soak and squeeze the sponge to activate D-Zyme™ Enzymatic Concentrate</p> <p>Options:</p> <ul style="list-style-type: none"> D-Zyme™ Concentrate Flat w/ 500 ml clear Basin - 100 ct / 1 case D-Zyme™ Concentrate Tubular w/ 500 ml clear Basin - 50 ct / 1 case <p><i>Also available:</i></p> <p>**D-Zyme™ Concentrate Flat in Master Bag & 500ml Basin - Add Water - 50 ct / 1 case</p> <p>**D-Zyme™ Concentrate Tubular in Master Bag & 500ml Basin - Add Water - 24 ct / 1case</p> | <p>MED-ETF500100</p> <p>MED-ETT500050</p> |  <p>D-Zyme™ Flat Concentrate & 500 ml Basin</p>  <p>D-Zyme™ Tubular Concentrate & 500 ml Basin</p> |
|--|---|--|



info@falconfireinc.com
www.medtrica.com

FALCONFIRE TECHNOLOGIES INC.

80 Bradford St., Suite WW23
Barrie, ON L4N 6S7

Toll Free: 866.635.2550